

Increased blood flow  
generates a  
natural healing response

geko™  
wound therapy

# Wound Therapy

Providing increased blood circulation to promote  
wound healing naturally from the inside



The geko™ Wound Therapy device **stimulates** the common peroneal nerve **activating** the calf and foot muscle pumps, increasing venous, arterial and microcirculatory blood flow.<sup>1, 2, 3</sup>

The blood flow increase prevents and reduces edema<sup>4</sup> and maintains TCpO<sub>2</sub> – promoting conditions favorable for wound healing.<sup>3, 5</sup>

For further information on the geko™ Wound Therapy device and the clinical evidence, please visit our website: [www.gekowound.ca](http://www.gekowound.ca)

## Demonstrations / Ordering:

T: +1 888 244 5579

Email: [orders@perfusemedtec.com](mailto:orders@perfusemedtec.com)

1. Tucker AT et al. Int J Angiol. 2010 Spring;19 (1): e31-e37
2. Williams KJ et al. Phlebology. 2015 Jun; 30 (5): 365-72
3. Jawad H et al. 2014 Journal Vasc Surg. Vol 2: 160-65
4. Williams KJ et al. Poster. Vascular Society Annual Scientific Meeting, Glasgow November 2014
5. Warwick D et al. Int J Angiol. 2015 April; 34 (2): 158-65

Perfuse  
Medtec